



The Spider Screw® System

The Spider Screw's geometry is a result of careful design in every single detail. In fact, the Spider Screw has obtained two international patents since its inception, due to its innovative characteristics: the simultaneous presence of the external and internal rectangular slots and round internal slots.

The Spider Screw's design is extremely versatile, due to its small dimensions and unique design. It is easily placed in either the maxilla or mandible, even where access is limited and bone quality is less than ideal. Placement is simplified by the self-drilling feature found in the K1 Plus, K2 Plus, K1 and K2 Spider Screw systems. Spider Pin requires pre-drilling.

The Spider Screw has been developed to offer a number of versatile anchorage options capable of immediate loading. Immediate loading is possible because the Spider Screw is a non-osteointegrable implant and consequently force can be applied immediately after placement. The applied force can range from 50 to 300 grams depending on screw choice, bone quality, and the desired orthodontic movement.

The Spider Screw is an anchorage device that can be used during every phase of orthodontic treatment and is suitable for symmetric or asymmetric anchorage. The Spider Screw assists in the success of orthodontic treatment, both in adults and adolescents, by reducing treatment times without patient compliance.

The Spider Screw package includes three removable labels containing important information (device name, reference code, lot number, etc.) which is to be applied to the patient's record card for traceability. Our sterile packaging ensures the Spider Screws are ready to use whenever you need them, saving your office valuable processing time. The screws have been cleaned, decontaminated and sterilized via gamma radiation, simply use sterile handling protocol.



Anchorage



Orthodontic Head

- Designed to facilitate appliance (wires, springs, etc.) placement
- Features two intersecting .022" slots
- Under tie-wing area also functions as another .022" x .025" slot featuring two intersecting slots of .027" in diameter
- Chamfered inlets simplify insertion of wires or ligatures
- Small head dimension is designed for patient comfort



Transmucosal Portion

- Variable length allows for optimal adaptation to different intraoral mucosa thicknesses during bio-maintenance
- Short length for areas of thinly attached gingiva
- Long for areas with thick or freely movable tissues
- Polished with a special treatment to help avoid soft tissue irritation and make cleaning easier to accomplish



Intrabony Portion

- Thread shape has an asymmetrical profile for easy placement
- Ensures maximum stability and avoidance of bone stress
- Self-drilling and self-tapping makes pre-drilling before insertion unnecessary dependent upon bone structure
- Easy to place while reducing the risk of root damage



Clinical Cases

- Intrusion Posterior Areas
- Direct Anchorage Uprighting and Molar Intrusion
- Class II Correction
- Lower Molar Protraction Indirect Anchorage
- Molar Uprighting and Intrusion

Visit our website for detailed case studies and before/after photos

Spider Screw®

Temporary Anchorage Device System

Made From High Quality Materials

- Grade 5 titanium construction
- Nickel-free for sensitive patients
- One per pack



Self-Ligating K1 Plus Series

Conical: Self-Drilling Thread

The Spider Screw Self-Ligating TAD - K1 is self-drilling and self-tapping. Due to the design of the conical thread, drilling is eliminated in most areas of the mouth. In areas of high bone density, it may be necessary to utilize the 1.1 mm drill provided to penetrate the cortical plate. The locking head is color coded pink for easy identification.



Long Neck – 2.0 mm	
Length	Item #
6.5 mm	SXL-1506
8.0 mm	SXL-1508
10.0 mm	SXL-1510

3.9 mm diameter head
1.5 mm diameter body

Self-Ligating K2 Plus Series

Conical: Self-Drilling Thread

The Spider Screw Self-Ligating TAD - K2 is self-drilling and self-tapping. Due to the design of the conical thread, drilling is eliminated in most areas of the mouth. In areas of high bone density, it may be necessary to utilize the 1.2 mm drill provided to penetrate the cortical plate. The locking head is color coded green for easy identification.



Long Neck – 2.0 mm	
Length	Item #
6.0 mm	SXL-1906
7.0 mm	SXL-1907
9.0 mm	SXL-1909
11.0 mm	SXL-1911

3.9 mm diameter head
1.9 mm diameter body

K1 Series

Conical: Self-Drilling Thread

The Spider Screw - K1 is self-drilling and self-tapping. Due to the design of the conical thread, drilling is eliminated in most areas of the mouth. In areas of high bone density, it may be necessary to utilize the 1.1 mm drill provided to penetrate the cortical plate.



Short Neck – 1.0 mm	
Length	Item #
6.5 mm	SCR-1506
8.0 mm	SCR-1508
10.0 mm	SCR-1510

3.4 mm diameter head
1.5 mm diameter body



Long Neck – 2.0 mm	
Length	Item #
6.5 mm	SCL-1506
8.0 mm	SCL-1508
10.0 mm	SCL-1510

3.4 mm diameter head
1.5 mm diameter body

K2 Series

Conical: Self-Drilling Thread

The Spider Screw - K2 is self-drilling and self-tapping. Due to the design of the conical thread, drilling is eliminated in most areas of the mouth. In areas of high bone density, it may be necessary to utilize the 1.2 mm drill provided to penetrate the cortical plate.



Short Neck – 1.0 mm	
Length	Item #
6.0 mm	SCR-1906
7.0 mm	SCR-1907
9.0 mm	SCR-1909
11.0 mm	SCR-1911

3.4 mm diameter head
1.9 mm diameter body



Long Neck – 2.0 mm	
Length	Item #
6.0 mm	SCL-1906
7.0 mm	SCL-1907
9.0 mm	SCL-1909
11.0 mm	SCL-1911

3.4 mm diameter head
1.9 mm diameter body

Spider Pin™

Cylindrical: Pre-Drilling Thread

The Spider Pin is self-tapping and requires pre-drilling. It is ideal for areas where a reduced size head is required and the infrabony portion is narrow. The simple design increases patient comfort and makes it easy to attach coil springs or elastic chains.



Long Neck – 2.0 mm	
Length	Item #
8.0 mm	SCL-1308
10.0 mm	SCL-1310

2.6 mm diameter head
1.3 mm diameter body

**Individual Self-Ligating
K1 & K2 Plus Series Screws**
Packaged sterile


Individual K1 & K2 Screws
Packaged sterile

Individual Spider Pins
Packaged sterile

Short Neck: Reduced neck height for thin tissue (anterior and lateral areas) Long Neck: Oversize neck height for soft thick tissue (posterior and lateral areas)



Spider Screw[®] Self-Ligating TAD



The locking head is color coded for easy identification and product selection



Head slot fits wires up to .018" x .022"

Long collared neck

Asymmetrical threads for easy placement and stability

Self-drilling and self-tapping with no pre-drilling required

Anchorage



Insert wire and turn Spider Screw Self-Ligating TAD head 90° with Self-Ligating Driver to secure

© 2021 Ortho Technology, Inc. Spider Screw is a registered trademark of Ortho Technology. Spider Screw Self-Ligating TAD is patented in the European Union and United States.



Spider Screw Self-Ligating TAD Starter Kit

- 1x 1.1 mm x 8.0 mm drill - K1
- 1x 1.3 mm x 9.0 mm drill - K2
- 1x Pick-up driver shaft - SXL
- 1x Contra angle pick-up driver short - SXL
- 1x Self-ligating TAD key
- 1x Handle driver - SXL
- 1x Screw driver body
- 1x Organizer

Entire Kit Can Be Sterilized

- Autoclavable up to 273°F/134°C
- Heat Sterilizable up to 356°F/180°C
- Chemically up to 297°F/132°C

Spider Screw Self-Ligating Starter Kit
Item #: CSS-6008

Replacement Organizer
Item #: CSS-6000

1 Year Warranty on Instruments



Built-in millimeter scale ensures correct Spider Screw selection



Included in Starter Kit



8.0 mm Drill - K1
Item #: FSC-1108-S
1.1 mm x 8.0 mm



9.0 mm Drill - K2
Item #: FSC-1309
1.3 mm x 9.0 mm



Pick-Up Driver Shaft - SXL
Item #: DSP-5652S
52.0 mm



Contra Angle Pick-Up Driver Short - SXL
Item #: DPQ-3420
20.0 mm



Self-Ligating TAD Key
Item #: DXL-2820



Handle Driver - SXL
Item #: DSQ-3424



Screw Driver Body
Item #: DSX-1690N-S

Optional Items



Contra Angle Pick-Up Driver Long - SXL
Item #: DPQ-3425
25.0 mm



Great for Larger Hands or for a **Firmer Grip**
Optional Round Replacement End for Screw Driver Body DSX-1690N-S
Item #: RCX-2545



Wide and Round End Handles provide easier, more ergonomic grip for larger hands
Screw Driver Body Wide Handle
Item #: DSX-1690RC



Spider Screw K1, K2 Starter Kit

- 1x 1.1 mm x 8.0 mm drill - K1
- 1x 1.3 mm x 9.0 mm drill - K2
- 1x Pick-up driver shaft
- 1x Cross driver shaft
- 1x Contra angle pick-up driver
- 1x Handle driver
- 1x Screw driver body
- 1x Organizer

Entire Kit Can Be Sterilized

- Autoclavable up to 273°F/134°C
- Heat Sterilizable up to 356°F/180°C
- Chemiclavable up to 297°F/132°C

Spider Screw Starter Kit

Item #: CSS-4008

Replacement Organizer - K1, K2

Item #: CSS-4000

1 Year Warranty on Instruments



Built-in millimeter scale ensures correct Spider Screw selection



Included in Starter Kit



8.0 mm Drill - K1
Item #: FSC-1108-S
1.1 mm x 8.0 mm



9.0 mm Drill - K2
Item #: FSC-1309
1.3 mm x 9.0 mm



Pick-Up Driver Shaft
Item #: DSP-5052N-S
52.0 mm



Cross Driver Shaft
Item #: DSX-2852N-S
52.0 mm



Contra Angle Pick-Up Driver
Item #: DPQ-2820-S
20.0 mm



Handle Driver
Item #: DSQ-2824-S



Screw Driver Body
Item #: DSX-1690N-S

Optional Items



- Automatically releases when the calibrated torque is reached to avoid overloading
- Torque can be set from 5 to 20 Ncm
- Helps eliminate screw breakage

Torque Screw Driver Body

Item #: DST-1600-S[«]

Note: Pickup driver sold separately. Sterilizable, see instructions for details.



10.0 mm Drill
Item #: FSC-1210-S[«]
1.2 mm x 10.0 mm



Contra Angle Pick-Up Driver
Item #: DPQ-2825
25.0 mm



Pick-Up Handle Driver
Item #: DPH-2824-S



Contra Angle Cross Driver
Item #: DPX-2830-S
30.0 mm



[«]Available by special order. Additional shipping time may be required.



Spider Pin™

Temporary Anchorage Device System

For Areas Where a Reduced Head Size is Required

- Simple and easy to use
- Smooth, rounded design for patient comfort
- Perfect for NiTi closed coil spring attachments
- Infrabony portion diameter is 1.3 mm at the widest point
- Ideal for narrow interproximal spaces
- Requires no patient cooperation and reduces treatment time
- Self-tapping and requires pre-drilling



Spider Pin Starter Kit

1.3 mm Diameter

- 1x 0.9 mm x 10.0 mm drill
- 1x Contra angle pick-up driver - pin
- 1x Short handle driver - pin
- 1x Pick-up handle driver shaft - pin
- 1x Screw driver body
- 1x Organizer

Entire Kit Can Be Sterilized

- Autoclavable up to 273°F/134°C
- Heat Sterilizable up to 356°F/180°C
- Chemoclavable up to 297°F/132°C

Spider Pin Starter Kit

Item #: CSS-3006

1 Year Warranty on Instruments



Included in Starter Kit



10.0 mm Drill - Pin
Item #: FSC-0910-S«
0.9 mm x 10.0 mm



Contra Angle Pick-Up Driver - Pin
Item #: DPQ-2322-S«
20.0 mm



Pick-Up Handle Driver Shaft - Pin
Item #: DSP-2352N-S«
52.0 mm



Short Handle Driver - Pin
Item #: DSQ-2324-S«



Screw Driver Body
Item #: DSX-1690N-S«
99.0 mm

© 2021 Ortho Technology, Inc. Spider Pin is a trademark of Ortho Technology. *Available by special order. Additional shipping time may be required.