

SPECIALS

BUR & POLISHER OFFER!

Buy 3 Get 1 FREE!*
Mix & Match!



High Speed Friction Grip Shank
1.6 mm width



Slow Speed Latch Type Shank
2.35 mm width



I Medium
J Fine
K Coarse

Debonding Carbide Burs

- A 12 Fluted Long Flame**
Item #: H48L-21-012 (Slow Speed)
Item #: H48L-31-012 (High Speed)
8.0 mm head length
- B 12 Fluted Football**
Item #: H379-21-023 (Slow Speed)
Item #: H379-31-023 (High Speed)
4.4 mm head length, 5 per pack
- C 12 Fluted Long Taper**
Item #: H375R-21-016 (Slow Speed)
Item #: H375R-31-016 (High Speed)
8.0 mm head length
- Slow Speed Latch Shank**
1 pack\$60.20/pk
3 or more\$56.20/pk
5 per pack
- High Speed Friction Grip Shank**
1 pack\$54.00/pk
3 or more\$51.00/pk
5 per pack
- D Terminator Metal Cutter**
Item #: T-2
1 pack\$28.20/pk
3 or more\$26.20/pk
High speed, 4.0 mm head length, 5 per pack

Dry heat sterilizable up to 340°F/170°C or autoclavable up to 250°F/121°C

Laboratory Carbide Burs

- E Medium Taper**
Item #: 79E-11-040
1 bur\$46.80 ea.
3 or more\$44.80 ea.
Low speed, 13.0 mm head length
- F Large Pointed Taper**
Item #: 257RE-11-060
1 bur\$46.80 ea.
3 or more\$44.80 ea.
Low speed, 14.0 mm head length
- G Large Egg Shape**
Item #: 251EF-11-060
1 bur\$46.80 ea.
3 or more\$44.80 ea.
Low speed, 14.0 mm head length
- H Vacuum Form Cutter Spiral Flute**
Item #: H219S-11-023
1 bur\$39.00 ea.
3 or more\$37.00 ea.
Low speed, 16.0 mm head length

Dry heat sterilizable up to 340°F/170°C or autoclavable up to 250°F/121°C



Polishing Wheels

- I 1670092 J 1670091 K 1670016
1 pack\$40.90/pk
5 or more\$36.90/pk
12 per pack

Finishing Polishing Points

- J High Shine Enamel Polisher 2.35 mm width**
Item #: W16D
1 pack\$45.10/pk
3 or more\$42.10/pk
Slow speed, Max. RPM 4,000-8,000. 10.4 mm head length, 3 per pack
- K High Shine Enamel Polisher 1.6 mm width**
Item #: W16F
1 pack\$45.10/pk
3 or more\$42.10/pk
High speed, Max. RPM 10,000. 8.0 mm head length, 3 per pack

© 2018 Ortho Technology, Inc. Galaxy is a trademark of Ortho Technology. *Buy 3 Galaxy products from this page get 1 of equal or lesser value free. Mix and match.

Galaxy™

INTERPROXIMAL REDUCTION SUPPLIES

BUY 3 GET 1 FREE!

Mix & Match! Buy any 3 Interproximal Reduction Supplies Get 1 of equal or lesser value FREE!



Diamond Disc Safety Guards

- A Metal Safety Guard**
19.0 mm 88611
22.0 mm 88612
1 disc guard\$70.20 ea.
3 or more\$66.20 ea.
- B Clear Safety Guard**
Item #: 88615
1 pack\$54.00/pk
3 or more\$51.00/pk
15 per pack, disposable

Galaxy IPR Diamond Strips

- Interproximal Strips Narrow**
..... 88603 88602 88601
Strip width 2.50 mm, 10 per pack
- Interproximal Strips Wide**
..... 88503 88502 88501
Strip width 3.75 mm, 10 per pack
- 1 pack\$114.20/pk
3 or more\$109.20/pk
- Interproximal Kits**
Narrow - 2.50 mm 88600
Wide - 3.75 mm 88500
1 kit\$107.20/kit
3 or more\$100.20/kit
3x extra fine, 3x fine, 3x medium

Diamond Disc Intro Kits

- Diamond 4 Disc - 22.0 mm Intro Kit**
Item #: DDKIT1
1 kit\$140.60 ea.
1x double-sided, 1x single-sided, 1x double-sided see-thru,
1x double-sided mesh
- Diamond 3 Disc - 22.0 mm Intro Kit**
Item #: DD-KIT
1 kit\$108.80 ea.
1x double-sided, 1x single-sided, 1x double-sided see-thru
- Diamond 3 Disc - 19.0 mm Intro Kit**
Item #: DDKIT2
1 kit\$108.80 ea.
1x double-sided, 1x single-sided, 1x double-sided mesh
All kits include mandrels and bur block

Galaxy Diamond Discs

- E Double-Sided Diamond Disc**
19.0 mm 88617
22.0 mm 88613
1 pack\$34.10/pk
3 or more\$31.10/pk
Disc thickness 0.19 mm
- F Double-Sided See-Thru Diamond Disc**
22.0 mm 88655
1 pack\$48.00/pk
3 or more\$46.00/pk
Disc thickness 0.20 mm
- G Double-Sided Mesh Diamond Disc**
19.0 mm 88657
22.0 mm 88656
1 pack\$48.00/pk
3 or more\$46.00/pk
Disc thickness 0.20 mm
- H Single-Sided Diamond Disc**
19.0 mm 88616
22.0 mm 88614
1 pack\$34.10/pk
3 or more\$31.10/pk
Disc thickness 0.15 mm, mounted diamond side down

© 2018 Ortho Technology, Inc. Galaxy is a trademark of Ortho Technology.