

Want to Simplify your Class II/III Treatment?

Treat Class II or III with reduced inventory and less breakage!

TruEase™

Class II/III Bite Corrector Devices



BUY 5 GET 2 FREE!*

• More Comfortable than Open-Spring Correctors

The fully enclosed spring is more hygienic and has no exposed mechanisms to impinge soft tissue

• No Reactivation Needed

"Out-of-the-box" installation in one appointment, and continuous forces until removal

• Simpler Inventory

Only two sizes to choose from per style - the universal design fits both sides of the mouth

• Also Treats Class III

Both Double Lock and Anchor Wire appliances can attach to either arch for Class II or III treatment

Class II Case Study

Patient with Class II, 1st division malocclusion and increased overbite and overjet.



No Bands & Headgear Tubes Needed



"The TruEase Double Lock version attaches directly onto orthodontic archwires without the necessity of laboratory work, and allows lateral mandibular movements."

- Dr. Pablo Echarri, Orthodontist



Bite Corrector Anchor Wire Kit

Standard Size (32.0-48.0 mm)..... TEBC210
Small Size (27.0-36.0 mm)..... TEBC211
1 pack **\$213.37/pk**
5 or more **\$204.10/pk**
2x appliances per kit

AS LOW AS **\$72.89**[†] each

Bite Corrector Double Lock Kit

Standard Size (27.0-36.0 mm)..... TEBC215
Small Size (23.0-32.0 mm)..... TEBC216
1 pack **\$239.19/pk**
5 or more **\$228.80/pk**
2x appliances per kit

AS LOW AS **\$81.71**[†] each



Call Us Today:
1-800-999-3161

© 2018 Ortho Technology, Inc. TruEase is a trademark of Ortho Technology. *Buy 5 packs of TruEase Bite Correctors get 2 packs free. †As low as price is based upon lowest advertised price including free product.


Exciting New Products to make Your Practice more Efficient & Profitable!

Buy 5 Packs Get 2 Free! Mix & Match

ALL ITEMS ON THIS PAGE

GALAXY™ DEBONDING BUR - 7901

New!



▪ Needle shaped trimming & finishing bur
▪ 12 fluted
Item #: 7901OT
4.0 mm head length, High Speed
1 pack **\$54.00/pk**
3 or more **\$51.00/pk**
5 per pack, Dry heat sterilizable up to 340°F/170°C or Autoclavable up to 250°F/121°C

AS LOW AS **\$729**[†] each

GALAXY™ IPR GAUGE SET

New!



Use during interproximal reduction to check space.
▪ Six Gauges Included:
.10 mm, .20 mm, .25 mm, .30 mm, .40 mm, & .50 mm
▪ Stainless steel
Item#: IPR-GAUGE
Autoclavable up to 248°F/120°C
1 pack **\$34.20 ea.**
3 or more **\$31.90 ea.**
1 per pack

AS LOW AS **\$2279**[†] each

ORTHOFLEX™ ROUND COMPOSITE BUTTONS

New!



Item #: OT-3300
2.9 mm one piece base
1 pack **\$50.20/pk**
3 or more **\$47.20/pk**
10 per pack

AS LOW AS **\$337**[†] each

RETAINER RETRIEVERS



1 pack **\$14.50/pk**
5 or more **\$13.80/pk**
10x patient individual-packs per pack

AS LOW AS **99¢**[†] each

Colors Available:
● ● ● ● ●

QWIKSTRIP™ SERRATED & ABRASIVE STRIPS



Single-Sided Assorted QS-ASST
Double-Sided Assorted..... QSDS-ASST
Curved Assorted QSC-ASST
Ultimate Assorted QS-ULT-10
Serrated .05 mm QS-SER
1 pack **\$91.40/pk**
3 or more **\$87.40/pk**
Steam sterilizable up to 325°F/162°C, Refill packs available 10 per pack

AS LOW AS **\$624**[†] each

CHEWIES™ ALIGNER TRAY SEATERS



▪ Helps seat aligner trays properly
▪ Chewing increases blood circulation
▪ Scented and unscented
▪ Made of soft Non-Latex** material
▪ 1 3/8" x 7/16"
1 pack **\$12.50/pk**
5 or more **\$10.50/pk**
10x patient 2 packs per pack

AS LOW AS **\$750**[†] /pack

© 2018 Ortho Technology, Inc. OrthoFlex is a trademark of Ortho Technology. †As low as price is based upon lowest advertised price including free product. *Free product must be of equal or lesser value.