

Thinking About Trying Self-Ligation?

Lotus Plus DS makes transitioning to self-ligation easy because it's designed like a traditional bracket. It lets you explore the efficiencies of an SLB, while allowing you to revert to elastomeric ligation and traditional mechanics if desired.

Lotus Plus DS

Self-Ligating Bracket System



Familiar Mini-Twin

Finish more accurately with easier placement and the option of ligating. Features a deep v-channel and clear vertical scribe line

Maximum Rotational Control

Doors run entire width of the archwire slot

Unique Planar, Thermal NiTi Clip

The passive & interactive clip designs reduce friction with flexible tabs, and withstand repeated opening for reliable forces throughout treatment

AS LOW AS
\$631[†]
/bracket



PASSIVE



HYBRID



INTERACTIVE

Doctors are making the switch from Damon® to finish cases on time. The NiTi door runs the entire width of the slot for maximum rotational control, and its flexible tabs create less friction for headache-free treatment.

Hybrid kits feature interactive brackets in the anterior and passive brackets in the posterior for enhanced mechanics.

In-Ovation® and Empower® users are switching to Lotus Plus DS to save time in their practice. The unique NiTi clip offers improved durability over standard chromium cobalt C-Clips and helps prevent doctors from having to re-bond brackets.

World's First Planar Interactive NiTi Clip!



With fewer stress points and a more resilient material, it effectively engages the archwire and does not deform as easily as chromium cobalt c-clips.

Available in Ortho Technology's version of the following prescriptions*:

Roth ■ MBT ■ Damon Super ■ Damon Standard ■ Damon Low

*Prescriptions are not implied to be an exact version of any other system.



Buy 4 Packs Get 1 Pack Free!

© 2019 Ortho Technology, Inc. Lotus Plus is a registered trademark of Ortho Technology. Domestic and international patents pending. Damon is a registered trademark of Ormco Corporation, In-Ovation is a registered trademark of GAC Dentsply, and Empower is a registered trademark of American Orthodontics.

World's First Thermal NiTi Clip Self-Ligating Buccal Tube!

Lotus Plus DS Self-Ligating Bracket System

14 Months • 4 Archwires • 6 Appointments



Beginning of Treatment



End of Treatment: 14 Months

Complete Control

"Interactivity is a key factor in SLB systems' performance. With Lotus Plus DS, you select the degree of engagement between the bracket and the archwire."

- Dr. Alex Lopes, Orthodontist

New! INTRODUCED AT THE AAO!

Self-Ligating Buccal Tubes

FOR A COMPLETE SL SYSTEM

Avoid step bends, decrease wire change time, and eliminate having to retie previously converted tubes with these innovative passive / interactive clip designs.

The bondable Lotus Plus DS Self-Ligating Buccal Tubes feature anatomical pads for maximum bond strength, and are a complete system for easier treatment.



INTRODUCTORY OFFER

Buy 3 Packs Get 1 Pack Free!

Lotus Plus DS Self-Ligating Buccal Tubes

1 pack \$126.30/pk
5 or more \$116.30/pk
1 pack (10 tubes)

AS LOW AS
\$872[†]
each

Available in Ortho Technology's version of the following prescriptions*:

Roth ■ Damon

*Prescriptions are not implied to be an exact version of any other system.



© 2019 Ortho Technology, Inc. †As low as price is based upon lowest advertised price including free product.