

# TIRED of being overcharged for METAL BRACKETS?

You might be paying for other companies' sales and marketing budgets.

The Competitor Price Tags = Huge Marketing Efforts and a large outside sales force with company cars, samples, and commissions.



AS LOW AS  
**\$201<sup>†</sup>**  
/bracket

## A Quality Bracket

Designed just a few years ago, Votion has some of the most up-to-date features for peace-of-mind treatment.

### COMPARE TO

OmniArch® PLUS by Dentsply Sirona

## Two-Piece Bracket System

### Designed for Clinical Efficiency and Patient Comfort

- Low profile design
- Ample under tie-wing area
- Anti-rotation available\*
- 80 gauge micro-etched mesh bonding base
- Intelligent hook design

### Available in the Most Popular Prescriptions

#### Single Patient Kit Description

	Slot Sizes
Upper/Lower 5x5	.018" .022"
Upper/Lower 5x5 Hooks on Cuspids	.018" .022"
Upper/Lower 5x5 Hooks 3, 4, & 5	.018" .022"
Upper/Lower 5x5 No Hooks Offset (U/L 4,5)	.018" .022"
Upper/Lower 5x5 Hooks on 3 Offsets (U/L 4,5)	.018" .022"
Upper/Lower 5x5 Hooks on 3, 4, 5 Offsets (U/L 4,5)	.018" .022"

#### Also Available: Individual Brackets

10 per pack



**Buy 20 Kits Get 10 Kits Free!**

**Call Us Today: 1-800-999-3161**

OmniArch is a registered trademark of Dentsply. \*Anti-Rotation available in Ortho Technology's version of the Roth prescription only.

**Buy 30 Get 10 Free!\***

or **Buy 100 Get 60 Free!\*** Mix and Match



Designed to Work Seamlessly with Votion® and Pinnacle® Bracket Systems

### Features:

- Lowered Hook
- Occlusal/Gingival Positioning Indents
- ID Color-Coded
- Metal Injection Molding (MIM) Stainless Steel
- Funneled Entrance  
With notch assists in easy archwire insertion into the buccal tube
- 80-Gauge Micro-Etched Mesh Pad  
Ensures maximum bond strength

Available in the Most Popular Prescriptions




AS LOW AS  
**\$424<sup>†</sup>**  
each

Defined Buccal Indent

### Votion® Bondable Buccal Tubes

Single Non-Convertible

#### Extended Bonding Base



Available Slot Sizes: .018" & .022"

1 pack (10 tubes)..... \$70.80/pk (\$7.08/tube)  
5 + packs (50 tubes)..... \$67.80/pk (\$6.78/tube)


Specifically Designed for the Most Comfortable Treatment



AS LOW AS  
**\$424<sup>†</sup>**  
each

### Votion® Mini Bondable Buccal Tubes

#### Mini Bonding Base



Available Slot Sizes: .018" & .022"

1 pack (10 tubes)..... \$70.80/pk (\$7.08/tube)  
5 + packs (50 tubes)..... \$67.80/pk (\$6.78/tube)

\*Buy 30 Votion Bondable Buccal Tubes or Mini Bondable Buccal Tubes get 10 free or buy 100 get 60 free. Mix and match. †As low as price is based upon lowest advertised price including free product.