

Tooth Tone® Archwires & Accessories

**BUY 5
GET 2
FREE!***

Mix & Match Wires



Stainless Steel Full Form Archwire

- Round wire sizes .014", .016", .018" and .020"
- Rectangle wire sizes .016" x .022", .017" x .025", .018" x .018", .019" x .025", .021" x .025", .020" x .020"

● Round	■ Rectangle
1 pack \$38.80/pk	1 pack \$39.70/pk
5 or more \$36.80/pk	5 or more \$37.70/pk
10 per pack	10 per pack

Nickel Titanium Full Form Archwire

- Round wire sizes .014", .016", .018" and .020"
- Rectangle wire sizes .016" x .022", .018" x .018", .019" x .025", .021" x .025", .020" x .020"

● Round	■ Rectangle
1 pack \$61.40/pk	1 pack \$65.90/pk
5 or more \$58.40/pk	5 or more \$62.90/pk
10 per pack	10 per pack

Note: Size shown is finished size. Tooth Tone coating adds .002" thickness to diameter. *Buy 5 packs of Tooth Tone Stainless Steel Full Form archwires or Tooth Tone Nickel Titanium Full Form archwires get 2 packs free. Mix and match. Free product must be of equal or lesser value.

- Superior Aesthetics
- Durable Coating
- Resistant to Staining
- Low Friction


Perfect With These Popular Aesthetic Brackets!



PURE®



Reflections™



Encore!™

Buy 4 Packs of Buttons Get 1 Pack of FREE!**



OrthoFlex™ Round Composite Buttons

Item #: OT-3300	1 pack \$50.20/pk
2.9 mm one piece base	3 or more \$47.20/pk
10 per pack	



OrthoFlex™ Rectangular Composite Buttons

Item #: OT-2344	1 pack \$50.20/pk
2.85 mm x 2.05 mm one piece base	3 or more \$47.20/pk
10 per pack	

Buy 4 Packs or Spools Get 1 Pack or Spool FREE!***

Make Aesthetic Treatment Complete With Elastomerics Non-Latex† MADE IN THE USA

Maximum™ Power Chain

	Closed	Short	Long
□ Clear	430-102	440-102	450-102
■ Tooth	430-103	440-103	450-103
□ Pearl	430-171	440-171	450-171
■ Blue Pearl	430-172	440-172	450-172

1 spool	\$26.20 ea.
10 spools.....	\$22.20 ea.
20 or more	\$18.20 ea.
<small>15' per spool (4.572 m)</small>	

Power Sticks™

Each Power Stick holds 50 molded ligatures. 1,000 per pack (20 sticks of 50), .115" ligatures

	Power Stick	Mini Power Stick
□ Clear (Stain Resistant)	404-002	407-002
■ Tooth	404-003	407-003
□ Pearl	404-071	407-071
■ Blue Pearl	404-072	407-072



Clear Elastomeric Power Stick

Pearl Elastomeric Mini Power Stick

Mini Power Sticks™

Each Mini Power Stick holds 10 molded ligatures to accommodate a single arch. Designed to reduce the potential risk of cross-contamination. 1,000 per pack (100 mini sticks of 10), .115" ligatures

1 pack	\$26.20/pk
10 packs.....	\$22.20/pk
20 or more	\$18.20/pk

© 2019 Ortho Technology, Inc. Tooth Tone and PURE are registered trademarks and OrthoFlex, Encore! and Reflections are trademark of Ortho Technology. *Buy 5 packs of Tooth Tone NiTi archwires get 2 packs of same item free. No mix and match. **Buy 4 packs of composite buttons get 1 of equal or lesser value free. ***Buy 4 packs or spools of Power Sticks or Power Chains get 1 of equal or lesser value free. Mix and match. †Not made from natural rubber latex.



THE SCIENCE OF PRECISION BONDING

Buy 3 Kits Get 1 Kit FREE!

Mix & Match! Buy any 3 Resilience Adhesive Kits from this page Get 1 of equal or lesser value Free!

Resilience®
Orthodontic Bonding Solutions



Light-Cure Band Cement

A light-activated band cement in a single paste system which provides a longer working time and a rapid setting time. Cement is easy to see due to the blue color. Contains fluoride to reduce the chance of decalcification.

AS LOW AS \$106.73[†] ea.

Light-Cure Band Cement Kit

Item #: FC-RESI-038	1 kit.....	\$151.30 ea.
	3 or more	\$142.30 ea.
	6x 5 gm push syringes	

Light-Cure Orthodontic Adhesive

The light-activated system provides the optimum properties for bonding ceramic, metal or plastic brackets to teeth. Delivers ideal handling characteristics which eliminate bracket flotation and allow for easy flash removal. Produces a quick setting, strong mechanical bond which debonds easily. Contains fluoride to help protect against potential decalcification.

Light-Activated Syringe Kit

Item #: FC-RESI-006	1 kit.....	\$135.60 ea.
	3 or more	\$127.60 ea.
	4x 3.5 gm syringes, 1x 7 mL sealant resin, 1x 4 gm green gel etchant syringe with 5x 22 gauge syringe tips, 2x 10 cavity mixing wells, 1x 50 sheet mixing pad (1 5/8" x 2 1/8"), 1x brush holder with 40x 3/4" brushes, and 10x 2" brushes	



AS LOW AS \$95.70[†] ea.



Low Viscosity Light-Cure Flowable Composite

Dependable, strong, and more resistant to shock and stress. Resilience Low Viscosity Flowable Composite can be used in many ways such as bonding lingual retainers, Invisalign® attachments, occlusal build-ups for opening a posterior bite, indirect bonding and much more!

LV Flowable Composite Kit

Item #: FC-RESI-007	1 kit.....	\$99.40 ea.
	3 or more	\$93.40 ea.
	1x 1.6 gm syringe green etchant, 3x 2 gm flowable composite syringes, 1x 3 mL light cure enamel bonding sealant resin, 25x 22 gauge syringe tips, 10x 2 1/2" bonding brushes, instructions	

© 2019 Ortho Technology, Inc. Resilience is a registered trademark of Ortho Technology. †As low as pricing is based on lowest advertised price including free product.