


Designed to Work Seamlessly with Votion® & Pinnacle® Bracket Systems



Buy 30 Get 10 Free!*

or **Buy 100 Get 60 Free!*** Mix and Match

VOTION® Bondable Buccal Tubes

AS LOW AS
\$4.24†
each



Votion® Bondable Buccal Tubes

Defined Buccal Indent



Votion® Mini Bondable Buccal Tubes

FEATURES:

- Lowered Hook
- Occlusal/Gingival Positioning Indents
- ID Color-Coded
- Metal Injection Molding (MIM) Stainless Steel
- Funneled Entrance
With notch assists in easy archwire insertion into the buccal tube
- 80-Gauge Micro-Etched Mesh Pad
Ensures maximum bond strength

Votion® Bondable Buccal Tubes

Single Non-Convertible

Quadrant	Slot Torque	Distal Offset	Color Code	BONDABLE - Extended Base	
				.018" Slot Size	.022" Slot Size
R Upper Right	-14°	14°	0	351-11-1620	351-12-1620
R Upper Left	-14°	14°	0	351-11-2620	351-12-2620
R Lower Right	-25°	4°	1	351-11-3620	351-12-3620
R Lower Left	-25°	4°	1	351-11-4620	351-12-4620
M Upper Right	-14°	10°	0	351-21-1620	351-22-1620
M Upper Left	-14°	10°	0	351-21-2620	351-22-2620
M Lower Right	-20°	0°	0	351-21-3620	351-22-3620
M Lower Left	-20°	0°	0	351-21-4620	351-22-4620
1 pack (10 tubes)				*\$70.80/pk (*\$7.80/tube)	
● 5+ packs (50 tubes)				*\$67.80/pk (*\$6.78/tube)	

R The Ortho Technology version of the Roth Rx**

M The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx**

Votion® Mini Bondable Buccal Tubes

Specifically Designed for the Most Comfortable Treatment

Quadrant	Slot Torque	Distal Offset	Color Code	Slot Size	
				.018"	.022"
R Upper Right	-10°	0°	0	352-11-1620	352-12-1620
R Upper Left	-10°	0°	0	352-11-2620	352-12-2620
R Lower Right	-25°	0°	1	352-11-3620	352-12-3620
R Lower Left	-25°	0°	1	352-11-4620	352-12-4620
M Upper Right	-14°	10°	0	352-21-1620	352-22-1620
M Upper Left	-14°	10°	0	352-21-2620	352-22-2620
M Lower Right	-10°	0°	0	352-21-3620	352-22-3620
M Lower Left	-10°	0°	0	352-21-4620	352-22-4620
1 pack (10 tubes)				*\$70.80/pk (*\$7.08/tube)	
● 5+ packs (50 tubes)				*\$67.80/pk (*\$6.78/tube)	

R The Ortho Technology version of the Roth Rx**

M The Ortho Technology version of the McLaughlin, Bennett, & Trevisi Rx**

**The Ortho Technology versions of the Roth and McLaughlin, Bennett, & Trevisi prescriptions are not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

© 2019 Ortho Technology, Inc. Votion is a registered trademark of Ortho Technology. *Buy 30 Votion Bondable Buccal Tubes or Mini Bondable Buccal Tubes get 10 free or buy 100 get 60 free. Mix and match. †As low as price is based upon lowest advertised price including free product.

Tired of WAITING WEEKS for your Band Assembly Order?

We have over 10,000 of the most popular Band and Tube Rx's ready for same day shipping.

TruFit™ 2.0

Prewelded Band and Tube Assemblies

OVER 10,000 IN-STOCK AND READY TO SHIP!

Place a first time order or restocking order of 250 or more TruFit™ 2.0 prewelded band and tube assemblies and **PRICE-LOCK™** at \$5.19 for all of 2019, including reorders of any quantity. **Guaranteed** through 12/31/19.

• Roth Rx* • McLaughlin, Bennett, & Trevisi Rx* • Standard Edgewise Rx



ONLY
\$5.19
Per Assembly

PRICE LOCK!

MADE IN THE USA

SWITCHING IS EASY With Our Band Size Conversion Trays



FREE!*

See our Band Conversion Chart Pages 66-67 in our Product Catalog

FEATURES AND BENEFITS:

- Precise, Tight Fit on Every Tooth
- Micro-Etched Interior Surface
- Laser Etched Exterior
- Sizes Available 29.5 up to 44
- Available With or Without Cleats

CUSTOM ORDERS AVAILABLE

Not Available for Same Day Shipment



Prewelded Sheath with Hook†
Special Offer Price \$2.00 ea. additional.
Prewelded Sheath†
Special Offer Price \$2.00 ea. additional.
Prewelded Button†
Special Offer Price \$1.00 ea. additional.



LOCK in Your Price for 2019!

Place a first time order or restocking order of 250 or more TruFit™ 2.0 prewelded band and tube assemblies and **PRICE-LOCK™** at \$5.19 for all of 2019, including reorders of any quantity. **Guaranteed** through 12/31/19.

**Call Us Today:
1-800-999-3161**

*The Ortho Technology versions of the Roth and McLaughlin, Bennett, & Trevisi prescriptions are not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

© 2019 Ortho Technology, Inc. TruFit is a trademark of Ortho Technology. †Lingual Sheaths are \$2.00 and lingual buttons are \$1.00 additional cost. Prices subject to change without notice. All bands are supplied with standard weld position unless custom ordered. Custom orders are final and cannot be returned. Allow 3-4 weeks for delivery. **Purchase 250 TruFit 2.0 Prewelded Band & Tube Assemblies, get 4 Stackable Band Organizers free.



Find us on **Facebook** | Follow us on **Twitter** | View us on **YouTube** | 11