



PURE[®]

Sapphire Bracket System



COMPARE TO
INSPIRE ICE™ by ORMCO

RISK FREE TRIAL OFFER!

Buy 5 Kits Get 2 Kits Free!

BUY 5 SINGLE PATIENT KITS (100 brackets)
GET 2 SINGLE PATIENT KITS (40 brackets) FREE!

MONEY BACK SATISFACTION GUARANTEED!

AS LOW AS
\$585[†]
ea.

Single Patient Kit Description	
Upper/Lower 5x5 1 kit (20 brackets)	\$183.80 kit (\$9.19 ea.)
● 10 kits (200 brackets)	\$163.80 kit (\$8.19 ea.)

ALSO AVAILABLE

Individual Brackets 10 brackets per pack
Upper 5-5 Single Patient Kit 10 brackets per kit
Upper 3x3 Single Patient Kit 12 brackets per kit

© 2020 Ortho Technology, Inc. PURE is a registered trademark of Ortho Technology.
†As low as price is based upon lowest advertised price including free product. Inspire ICE is a trademark of Ormco/Sybron Dental Specialties.

THE ULTIMATE in Crystal Clear Aesthetics™

- Smooth, strong and contoured tie-wings
- Aluminum oxide spheres on contoured base
- Ultra-smooth, strong hooks
- Removable color coded ID dot
- Ample tie-wing area

AVAILABLE IN ORTHO TECHNOLOGY'S VERSION OF THESE PRESCRIPTIONS*


- Roth* • McLaughlin, Bennett, & Trevisi*
- Standard Edgewise

*The Ortho Technology versions of the Roth and McLaughlin, Bennett, & Trevisi prescriptions are not implied to be an exact version of any other system, nor do we claim any endorsement of the respective doctors.

Tooth Tone[®] Archwires & Accessories

BUY 5 GET 2 FREE!*

Mix & Match Wires



- Superior Aesthetics
- Durable Coating
- Resistant to Staining
- Low Friction

Perfect With These
Popular Aesthetic Brackets!



Nickel Titanium Full Form Archwire

● Round wire sizes .014", .016", .018" and .020"	
■ Rectangle wire sizes .016" x .022", .018" x .018", .019" x .025", .021" x .025", .020" x .020"	
● Round 1 pack \$63.20/pk 5 or more \$60.20/pk 10 per pack	■ Rectangle 1 pack \$67.80/pk 5 or more \$64.80/pk 10 per pack

Stainless Steel Full Form Archwire

● Round wire sizes .014", .016", .018" and .020"	
■ Rectangle wire sizes .016" x .022", .017" x .025", .018" x .018", .019" x .025", .021" x .025", .020" x .020"	
● Round 1 pack \$39.90/pk 5 or more \$37.90/pk 10 per pack	■ Rectangle 1 pack \$40.80/pk 5 or more \$38.80/pk 10 per pack

Note: Size shown is finished size. Tooth Tone coating adds .002" thickness to diameter.
*Buy 5 packs of Tooth Tone Stainless Steel Full Form archwires or Tooth Tone Nickel Titanium Full Form archwires get 2 packs free. Mix and match. Free product must be of equal or lesser value.

Buy 5 Packs of Buttons Get 1 Pack of FREE!**



OrthoFlex™ Round Composite Buttons	
Item #: OT-3300 2.9 mm one piece base 10 per pack	1 pack \$51.60/pk 3 or more \$48.60/pk



OrthoFlex™ Rectangular Composite Buttons	
Item #: OT-2344 2.85 mm x 2.05 mm one piece base 10 per pack	1 pack \$51.60/pk 3 or more \$48.60/pk

Buy 4 Packs or Spools Get 1 Pack or Spool FREE!***

Make Aesthetic Treatment Complete With Elastomerics

Non-Latex† **MADE IN THE USA**


Maximum™ Power Chain

	Closed	Short	Long
□ Clear	430-102	440-102	450-102
■ Tooth	430-103	440-103	450-103
□ Pearl	430-171	440-171	450-171
■ Blue Pearl	430-172	440-172	450-172
1 spool	\$26.70 ea.		
10 spools.....	\$22.70 ea.		
20 or more	\$18.70 ea.		
15' per spool (4.572 m)			

Power Sticks™

Each Power Stick holds 50 molded ligatures.
1,000 per pack (20 sticks of 50), .115" ligatures

	Power Stick	Mini Power Stick
□ Clear (Stain Resistant)	404-002	407-002
■ Tooth	404-003	407-003
□ Pearl	404-071	407-071
■ Blue Pearl	404-072	407-072



Mini Power Sticks™

Each Mini Power Stick holds 10 molded ligatures to accommodate a single arch. Designed to reduce the potential risk of cross-contamination.
1,000 per pack (100 mini sticks of 10), .115" ligatures

1 pack	\$26.70/pk
10 packs.....	\$22.70/pk
20 or more	\$18.70/pk

© 2020 Ortho Technology, Inc. Tooth Tone and PURE are registered trademarks and OrthoFlex, Encore! and Reflections are trademark of Ortho Technology. **Buy 5 packs of composite buttons get 1 pack free. ***Buy 4 packs or spools of Power Sticks or Power Chains get 1 free. Mix and match. †Not made from natural rubber latex.