

Thinking About Trying Self-Ligation?

Lotus Plus DS makes transitioning to self-ligation easy because it's designed like a traditional bracket. It lets you explore the efficiencies of an SLB, while allowing you to utilize elastomeric ligation and traditional mechanics if desired.

Lotus PLUS DS

Self-Ligating Bracket System

AS LOW AS
\$650†
/bracket

Familiar Mini-Twin

Finish more accurately with easier placement and the option of ligating. Features a deep v-channel and clear vertical scribe line

Maximum Rotational Control

Doors run entire width of the archwire slot

Unique Planar, Thermal NiTi Clip

The passive & interactive clip designs reduce friction with flexible tabs, and withstand repeated opening for reliable forces throughout treatment



Dual System, Unlimited Options

- Passive
- Interactive
- Hybrid



PASSIVE

Doctors are making the switch from Damon® to finish cases on time. The NiTi door runs the entire width of the slot for maximum rotational control, and its flexible tabs create less friction during treatment.



HYBRID

Hybrid kits feature interactive brackets in the anterior and passive brackets in the posterior for enhanced mechanics.



INTERACTIVE

In-Ovation® and Empower® users are switching to Lotus Plus DS to save time in their practice. The unique NiTi clip offers improved durability over standard chromium cobalt C-clips and helps prevent doctors from having to re-bond brackets.

World's First Planar Interactive NiTi Clip!



With fewer stress points and a more resilient material, it effectively engages the archwire and does not deform as easily as chromium cobalt c-clips.

Available in Ortho Technology's version of the following prescriptions*:

Roth • MBT • Damon Super • Damon Standard • Damon Low

*Prescriptions are not implied to be an exact version of any other system.

Buy 4 Kits Get 1 Kit Free!

© 2020 Ortho Technology, Inc. Lotus Plus is a registered trademark of Ortho Technology. Domestic and international patents pending. Damon is a registered trademark of Ormco Corporation, In-Ovation is a registered trademark of GAC Dentsply, and Empower is a registered trademark of American Orthodontics. †Buy 4 kits get 1 kit free. As low as price is based upon retail price including free product.

World's First Thermal NiTi Planar Clip Self-Ligating Buccal Tube!

Lotus PLUS DS Self-Ligating Bracket System

14 Months • 4 Archwires • 6 Appointments



Beginning of Treatment



End of Treatment: 14 Months

Complete Control

"Interactivity is a key factor in SLB systems' performance. With Lotus Plus DS, you select the degree of engagement between the bracket and the archwire."

- Dr. Alex Lopes, Orthodontist

Self-Ligating Buccal Tubes FOR A COMPLETE SL SYSTEM

An Advanced Buccal Tube

Avoid step bends, decrease wire change time, and eliminate having to retie previously converted tubes with these innovative passive / interactive clip designs.

The bondable Lotus Plus DS Self-Ligating Buccal Tubes feature anatomical pads for maximum bond strength, and are a complete system for easier treatment.



INTRODUCTORY OFFER

Buy 3 Packs Get 1 Pack Free!

Lotus Plus DS Self-Ligating Buccal Tubes

1 pack \$129.80/pk
5 or more \$119.80/pk
1 pack (10 tubes)

AS LOW AS
\$899†
ea.

Available in Ortho Technology's version of the following prescriptions*:

Roth • Damon

*Prescriptions are not implied to be an exact version of any other system.

© 2020 Ortho Technology, Inc. †As low as price is based upon lowest advertised price including free product.

