

In-Ovation R Brackets

Lotus+ DS Brackets

Roth (Base) RX .022 Slot			Roth RX .022 Slot			
GAC ITEM #	Torque Angulation	Tooth	Tooth	OT ITEM #	Torque Angulation	Difference
Maxillary			Maxillary			
189-112-00	+12T +5A	UR1	UR1	LPAR22-001	+12T +5A	
189-212-00	+12T +5A	UL1	UL1	LPAR22-002	+12T +5A	
189-122-00	+8T +9A	UR2	UR2	LPAR22-011	+8T +9A	
189-222-00	+8T +9A	UL2	UL2	LPAR22-012	+8T +9A	
189-132-10	-2T +13A	UR3/HK	UR3/HK	LPAR22-021	0T +13A	2T, offset
189-232-10	-2T +13A	UL3/HK	UL3/HK	LPAR22-022	0T +13A	2T, offset
189-132-00	-2T +13A	UR3	UR3	LPAR22-023	0T +13A	2T, offset
189-232-00	-2T +13A	UL3	UL3	LPAR22-024	0T +13A	2T, offset
189-142-00	-7T 0A	UR4/5	UR4/5	LPAR22-033	-7T 0A	offset
189-242-00	-7T 0A	UL4/5	UL4/5	LPAR22-034	-7T 0A	offset
189-142-90	-7T 0A	UR4/5 HK MESIAL	UR4/5 HK	LPAR22-031	-7T 0A	offset
189-242-90	-7T 0A	UL4/5 HK MESIAL	UL4-5 HK	LPAR22-032	-7T 0A	offset
189-142-10	-7T 0A	UR4/5HK DISTAL	UR4/5 HK	LPAR22-031	-7T 0A	Mesial hook only, offset
189-242-10	-7T 0A	UR4/5 HK DISTAL	UL4/5 HK	LPAR22-032	-7T 0A	Mesial hook only, offset
87-162-05	-10T 0A	UR6 HK SL	UR6 HK	301-12-1020	-10T 0A	
87-262-05	-10T 0A	UL6 HK SL	UL6 HK	301-12-2620	-10T 0A	
Mandibular			Mandibular			
189-312-00	-1T 2A	LR1-2	LR1-2	LPAR22-040	0T 0A	1T 2A
189-412-00	-1T 2A	LL1-2	LL1-2	LPAR22-040	0T 0A	1T 2A
189-332-10	-11T +7A	LR3/HK	LR3/HK	LPAR22-051	-11T +5A	2T, offset
189-432-10	-11T +7A	LL3/HK	LL3/HK	LPAR22-052	-11T +5A	2T, offset
189-332-00	-11T +7A	LR3	LR3	LPAR22-053	-11T +5A	2T, offset
189-432-00	-11T +7A	LL3	LL3	LPAR22-054	-11T +5A	2T, offset
189-342-00	-17T -1A	LR4	LR4	LPAR22-063	-17T 0A	1A, offset
189-442-00	-17T -1A	LL4	LL4	LPAR22-064	-17T 0A	1A, offset
87-342-10	-17T -1A	LR4/HK	LR4/HK	LPAR22-061	-17T 0A	1A, offset
87-442-10	-17T -1A	LL4/HK	LL4/HK	LPAR22-062	-17T 0A	1A, offset
189-342-08	-17T -1A	LR4/OB	LR4/Offset Pad	LPAR22-063	-17T 0A	1A
189-442-08	-17T -1A	LL4/OB	LL4/Offset Pad	LPAR22-064	-17T 0A	1A
189-352-00	-22T -1A	LR5	LR5	LPAR22-073	-22T 0A	1A, offset
189-452-00	-22T -1A	LL5	LL5	LPAR22-074	-22T 0A	1A, offset
189-352-10	-22T -1A	LR5/HK	LR5/HK	LPAR22-071	-22T 0A	1A, offset
189-452-10	-22T -1A	LL5/HK	LL5/HK	LPAR22-072	-22T 0A	1A, offset
189-352-08	-22T -1A	LR5/OB	LR5/Offset Pad	LPAR22-073	-22T 0A	1A
189-452-08	-22T -1A	LL5/OB	LL5/Offset Pad	LPAR22-074	-22T 0A	1A
87-362-09	-25T 0A	LR6 HK SL	LR6 HK	301-12-3620	-25T 0A	
87-462-09	-25T 0A	LL6 HK SL	LL6 HK	301-12-4620	-25T 0A	
KITS			KITS			
K189-055-23	n/a	U/L 5-5 Cuspid HK	U/L 5-5 Cuspid HK	LPAR22-1300	n/a	
K189-055-24	n/a	U/L 5-5 Cuspid Biscupid HKM	U/L 5-5 HKS 345	LPAR22-1345	n/a	
K189-055-22	n/a	U/L 5-5 Cuspid Biscupid HKD	UL 5-5 HKS 345	LPAR22-1345	n/a	Upper Cuspid hook distal, biscuspid hook mesial

ROTH 022

****Please note all cuspids and biscuspids on Lotus+ have offset pads****