

GAC OmniArch Plus Brackets

Ortho Technology Votion Brackets

ROTH RX 022

Roth RX .022 Slot			Roth RX .022 Slot				
GAC ITEM #	Torque Angulation	Tooth	Tooth	OT ITEM #	Torque Angulation	Difference	VOTION Anti-Rotation
Maxillary			Maxillary				
J77-112-00	+12T +5A	UR1	UR1	350-92-1100	+12T +5A		Centrals & Laterals = 0°
J77-212-00	+12T +5A	UL1	UL1	350-92-2100	+12T +5A		Centrals & Laterals = 0°
J77-122-00	+8T +9A	UR2	UR2	350-92-1200	+8T +9A		Centrals & Laterals = 0°
J77-222-00	+8T +9A	UL2	UL2	350-92-2200	+8T +9A		Centrals & Laterals = 0°
J77-132-10	+2T +10A	UR3/HK	UR3/HK	350-92-1311	-2T +13A	4T 3A	Cuspids = +4°M
J77-232-10	+2T +10A	UL3/HK	UL3/HK	350-92-2311	-2T +13A	4T 3A	Cuspids = +4°M
J77-542-00	-7T 0A	U4 Universal	UR4/5	350-92-1902	-7T 0A		Bicuspids = +2°D
			UL4/5	350-92-2902	-7T 0A		Bicuspids = +2°D
J77-142-10	-7T 0A	UR4/HK	UR4-5/HK	350-92-1912	-7T 0A		Bicuspids = +2°D
J77-242-10	-7T 0A	UL4/HK	UL4-5/HK	350-92-2912	-7T 0A		Bicuspids = +2°D
J77-542-07	-7T 0A	U4-5/OB (Offset Base)	UR4-5/Offset Pad	350-92-1922	-7T 0A		Bicuspids = +2°D
			UL4-5/Offset Pad	350-92-2922	-7T 0A		Bicuspids = +2°D
J77-552-00	-7T 0A	U5 Universal	UR4/5	350-92-1902	-7T 0A		Bicuspids = +2°D
			UL4/5	350-92-2902	-7T 0A		Bicuspids = +2°D
J77-152-10	-7T 0A	UR5/HK	UR4-5/HK	350-92-1912	-7T 0A		Bicuspids = +2°D
J77-252-10	-7T 0A	UL5/HK	UL4-5/HK	350-92-2912	-7T 0A		Bicuspids = +2°D
Mandibular			Mandibular				
J77-612-00	-1T 0A	L1-2	L1-2	350-92-6800	-1T 0A		Lower Anterior = 0°
J77-332-10	-11T +7A	LR3/HK	LR3/HK	350-92-3311	-11T +7A		Cuspids = +2°M
J77-432-10	-11T +7A	LL3/HK	LL3/HK	350-92-4311	-11T +7A		Cuspids = +2°M
J77-642-00	-17T 0A	L4 Universal	LR4	350-92-3402	-17T 0A		Bicuspids = +4°D
			LL4	350-92-4402	-17T 0A		Bicuspids = +4°D
J77-342-10	-17T 0A	LR4/HK	LR4/HK	350-92-3412	-17T 0A		Bicuspids = +4°D
J77-442-10	-17T 0A	LL4/HK	LL4/HK	350-92-4412	-17T 0A		Bicuspids = +4°D
J77-642-07	-17T 0A	L4/OB Universal	LR4/Offset Pad	350-92-3422	-17T 0A		Bicuspids = +4°D
			LL4/Offset Pad	350-92-4422	-17T 0A		Bicuspids = +4°D
J77-652-00	-22T 0A	L5 Universal	LR5	350-92-3502	-22T 0A		Bicuspids = +4°D
			LL5	350-92-4502	-22T 0A		Bicuspids = +4°D
J77-352-10	-22T 0A	LR5/HK	LR5/HK	350-92-3512	-22T 0A		Bicuspids = +4°D
J77-452-10	-22T 0A	LL5/HK	LL5/HK	350-92-4512	-22T 0A		Bicuspids = +4°D
J77-652-07	-22T 0A	L5/OB Universal	LR5/Offset Pad	350-92-3522	-22T 0A		Bicuspids = +4°D
			LL5/Offset Pad	350-92-4522	-22T 0A		Bicuspids = +4°D
KITS			KITS				
KJ77-055-23	n/a	U/L 5-5 Cuspid HKM	U/L 5-5 Cuspid HK (Distal)	350-92-2010	n/a		
KJ77-055-22	n/a	U/L 5-5 Cuspid BiCuspid HKM	U/L 5-5 Cuspid (HKD) BiCuspid (HKM)	350-92-2031	n/a		
TYPJ77-777-77	n/a	Consultation Model/Typodont	Consultation Model/Typodont	G50-122	n/a		