

Looking for a Modern Two-Piece Bracket?

Designed just a few years ago, Votion has some of the most up-to-date features for more efficient treatment, without the competitor price tag!

VOTION[®]
Stainless Steel Bracket System

The Bracket Built for Practice Productivity

COMPARE TO:
OmniArch[®] Plus by Dentsply Sirona



ORTHO TECHNOLOGY[®] Two-Piece Bracket System

Make the Switch from GAC & Get

20% OFF!*


The Votion Bracket System utilizes Metal Injection Molding (MIM) technology to achieve extremely complex bracket structures and contours that enhance treatment efficiency.

Intelligent Hook Design
Hooks are angled at 15° to reduce interference with gingival tissue

Smooth Facial Contours
For increased patient comfort

Vertical Scribe Line
Runs entire length of bracket base and pad for easy and more accurate placement on the long axis

Compound Contour Base
Superior fit and greater contact surface area for improved retention



Permanent Identification System
With a removable color coded dot and permanently recessed I.D. mark


Ample Under Tie-Wing Area
Provides easy and secure ligation

Rounded Archwire Slot Contours
Mesial and distal slot entry points are rounded to prevent archwire notching when rotations are present and during retraction/protraction procedures

Anti-Rotation

Ortho Technology's version of the Roth Prescription features built-in anti-rotation in the cuspid and bicuspid bracket archwire slots to prevent unwanted rotation during sliding mechanics

Call Today and Ask About
Our Bracket Exchange Program



Try Our Votion Bondable Buccal Tubes
Designed to Work Seamlessly with
Votion Bracket System!

© 2020 Ortho Technology, Inc. Votion is a registered trademark of Ortho Technology. OmniArch is a registered trademark of Dentsply GAC / Dentsply Sirona.
*Make the switch from Dentsply GAC products to Ortho Technology products, and receive 20% OFF your entire order.

IT'S EASY To Switch ARCHWIRES FROM GAC!

By offering an extensive selection of high quality wire types and archforms to treat every case, Ortho Technology is your one-stop shop!

SWITCH FROM GAC & GET 20% OFF!*



Ortho Technology Archwire Brands

Compare To: GAC Archwire Brands

TruFlex™ Thermal Nickel Titanium	Sentalloy® & Neo Sentalloy®
TruFlex™ Copper Nickel Titanium	Copperloy® Nickel Titanium
TruForce™ Stainless Steel	Stainless Steel
BetaForce™ Beta Titanium	Resolve® Beta Titanium
TRIFECTA® Thermal Nickel Titanium	BioForce® PLUS
Tooth Tone® Coated NiTi & Stainless Steel	Wires with IonGuard®
TruGold™ 24K Gold Plated	N/A



MADE IN THE USA

ORTHO TECHNOLOGY

CONVERSIONS Made SIMPLE!

We have thousands of the most popular Band and Tube Rx's ready for same day shipment.

TruFit™ 2.0

IN-STOCK & READY TO SHIP!

Prewelded Band and Tube Assemblies

MADE IN THE USA

Regardless of the brand currently being used in your practice, when you convert to TruFit 2.0, we will label our band trays with your current sizing, eliminating the need to learn our sizes and leaving your workflow unchanged.



COMPARE TO: Snap-Fit Molar Bands by GAC Dentsply

ORTHO TECHNOLOGY

Call Today and Ask About Our Band Exchange Program!

FREE Band Size Conversion Trays Included with Your Order



See our Band Conversion Chart Pages 66-67 in our Product Catalog

© 2020 Ortho Technology, Inc. TruFlex, TruForce, BetaForce, TruGold, and TruFit are trademarks and ToothTone and TRIFECTA are registered trademarks of Ortho Technology. Sentalloy, Copperloy, BioForce, and Resolve are registered trademarks of Dentsply GAC / Dentsply Sirona. *Make the switch from Dentsply GAC products to Ortho Technology products, and receive 20% OFF your entire order.

ORTHO TECHNOLOGY

US: 1.800.999.3161

Outside US: 1.813.501.1650

www.orthotechnology.com