



Forehead Rest with Ergonomic Dual Free-Flow Air Vents



Vertical Main Frame



Metal Reinforced Horizontal Crossbar



Durable Screws Throughout

Chin Cup with Ergonomic Streamlined Design, Free-Flow Air Vents and Full Range of Motion Option



Small Light-Weight End Caps

Multi-Adjustable Facemask®

Discover The Ideal Choice to Achieve Maximum Class III Compliance

- Free-flow air vents
- Full range of motion option
- Maximum patient comfort
- Fully adjustable
- Light-weight
- Screws designed for increased durability


MADE IN THE USA

Exclusive Patient Comfort Design

Quick and Easy Fit

1. Position bend in vertical bar over nose
2. Center forehead rest on forehead and lock in position
3. Center chin cup on chin and lock in position
4. Place crossbar approximately 2.0 mm below the upper lip and lock in position

KEY FEATURES



Free-Flow Air Vent

Forehead rest design features dual air vents for maximum comfort and air circulation.




Free-Flow Air Vent

Chin cup design features air vents for maximum comfort and air circulation.



Motion Option

Sliding design allows for a dynamic range of motion without side-to-side rotation.



American Pride
Item #: 550-406



Canadian Pride
Item #: 550-411



Green Camo
Item #: 550-413



Hawaiian Hibiscus
Item #: 550-407



Pink Camo
Item #: 550-412



Racin' Flames
Item #: 550-408



Union Jack
Item #: 550-410



Wild Cat Print
Item #: 550-409



Blue
Item #: 550-400



Clear
Item #: 550-403



Pearl White
Item #: 550-404



Red
Item #: 550-401



Teal
Item #: 550-405





Multi-Adjustable Facemask
Includes: 1x facemask, 1x hexagon key, and 2 sets of adhesive chin and forehead pads

See page 107 for Extraoral Elastics

© 2021 Ortho Technology, Inc.
Multi-Adjustable Facemask is a registered trademark of Ortho Technology.

Replacement Forehead Rest and Chin Cup Kit

American Pride	406-RK2
Canadian Pride.....	411-RK2
Green Camo.....	413-RK2
Hawaiian Hibiscus.....	407-RK2
Pink Camo	412-RK2
Racin' Flames.....	408-RK2
Union Jack	410-RK2
Wild Cat Print.....	409-RK2
Solid Colors	550-RK2
1x hexagon key and 2 sets of adhesive chin and forehead pads	



Large



Small

Reverse Pull Facemask

Fully Adjustable Comfort Design Aids in Patient Compliance

- Maximum patient comfort
- Treats maxillary insufficiencies
- Fully adjustable
- Stainless steel and plastic
- Leather padded contact points
- Hexagon key included
- Storage bag included

Large REV-PULL
 Small REV-PULS
 1x facemask, 1x hexagon key,
 1x storage bag
See page 105 for Extraoral Elastics

Replacement Pad Kit
 Item #: 100002
 1x forehead pad and 1x chin pad



A



B



Safety Release Modules

- Flexible tear-proof strap
- Durable steel clamp or easy to install plastic release clamp design
- Attaches to any neckpad
- Safety release
- Large facebow holes
- Royal blue colored

A. Steel C-Clamp
 600 gm Medium 550-110
 750 gm Heavy 550-111
 10 per pack

B. Plastic C-Clamp
 600 gm Medium 550-100
 750 gm Heavy 550-101
 10 per pack



Assorted Colors

Cervical Neckpads

- Comfortable cotton fabric
- Sturdy washable design
- Three module attachment locations
- Fun, stylish colors
- Use with any style releasable force module

Royal Blue..... 550-700
 Assorted Colors*..... 550-600
 5 per pack

*Royal Blue, Black, Denim, Crimson, Surf



Royal Blue



Royal Blue

High Pull Head Caps

- Durable nylon
- Sturdy washable design
- Two hook-up tabs
- Use with any style releasable force module

*Royal Blue, Navy Blue, Surf, Crimson, Black

Royal Blue..... 550-750
 Assorted Colors*..... 550-650
 5 per pack



Ortho Performance® Mesh Bags

- Headgear storage
- One bag with multiple uses
- Can be used to hold headgear items, multi-adjustable facemasks or oral hygiene kits
- Solid blue in color with a secure zipper closure


- Made from premium micro mesh material

Royal Blue.....560111
 10 per pack
 Measures 8 3/4" x 5"

© 2021 Ortho Technology, Inc. Ortho Performance is a registered trademark of Ortho Technology.



Full Form Extraoral Facebows



High Quality Material

High quality stainless steel is used in the manufacturing process of our facebows.

Wide Inner Bow Size Measurement

For accurate selection of the proper size, measurements are based on center of loop to center of opposite loop of the inner bow.

#1	83.0 mm	551-111
#2	90.0 mm	551-112
#3	97.0 mm	551-113
#4	104.0 mm	551-114
#5	111.0 mm	551-115

Universal (*inner bow without loops*)

Inner Bow Shapes

Long 551-110

Standard 551-210

1 per pack

TRY OUR...
TruFit™ 2.0
Molar Band
Assemblies
See pages 59-62

Short Outer Arch Facebows



- Durable brazed joint for unsurpassed strength
- Preformed omega loops on innerbow wire for easy adjustment

#1	83.0 mm	M4835-01
#2	90.0 mm	M4835-02
#3	97.0 mm	M4835-03
#4	104.0 mm	M4835-04
#5	111.0 mm	M4835-05

1 per pack