





**PixSure
Perfect**[®]

Clinical Photography Supplies

**Intraoral Rhodium Coated
Photographic Mirrors**

MADE
IN THE **USA**



**The Durable Surface Coating Produces
a Sharp and Precise Image Every Time**

A. Occlusal (Child)

Item #: 32700
3 1/3" (84.6 mm) x 2 2/5" (61.5 mm) x 2" (50.8 mm)

B. Occlusal (Child Extra Long)

Item #: 32701
4 1/2" (114.3 mm) x 2 1/2" (63.5 mm) x 2 1/2" (57.2 mm)

C. Occlusal (Adult)

Item #: 32702
4" (101.6 mm) x 3" (76.2 mm) x 2 1/2" (63.5 mm)

D. Occlusal (Adult Extra Large)

Item #: 32703
5 1/2" (139.7 mm) x 3" (76.2 mm) x 2 2/3" (67.8 mm)

E. Lingual

Item #: 32704
6" (152.4 mm) x 1 2/5" (35.6 mm) x 1 1/2" (38.1 mm)

F. Buccal (Narrow)

Item #: 32705
5 3/4" (145.0 mm) x 2" (50.8 mm) x 1" (25.4 mm)

G. Buccal

Item #: 32706
6 1/4" (158.8 mm) x 2 1/4" (57.2 mm) x 1 3/4" (44.4 mm)



(C) Occlusal (Adult)




(G) Buccal



(E) Lingual

- Superior color and contrast
- Rhodium coated glass
- Tempered glass coated on both sides
- Optimum reflectivity and adhesion
- Resistant to abrasion
- Extremely durable coating
- Use with digital or 35 mm
- Withstands repeated sterilization

1 per pack
Autoclavable up to 270°F/132°C



The surface of PixSure Perfect[®] rhodium coated mirrors ensures maximum reflectivity, adhesion, and resistance to abrasion

© 2021 Ortho Technology, Inc. PixSure Perfect is a registered trademark of Ortho Technology.





All PixSure Perfect® Intraoral Stainless Steel Double Ended Photographic Mirror surfaces are set at 30° to each other for ease of use



PixSure Perfect®

Clinical Photography Supplies

Intraoral Stainless Steel Photographic Mirrors

These Highly Polished Stainless Steel Double-Ended Mirrors Provide the Perfect Imaging Solution for Every Patient

- No double reflections
- Polished stainless steel
- Double ended versatility
- Extremely durable finish
- Use with digital or 35 mm
- Withstands repeated sterilization
- Beveled edge for patient comfort

The highly polished surfaces of PixSure Perfect® stainless steel mirrors ensure high reflectivity and are virtually indestructible

*1 per pack
Autoclavable up to 273°F/134°C*

- A. Adult Occlusal/Child Occlusal**.....32800
2³/₄" (70 mm) x 7" (178 mm) x 2¹/₄" (57 mm)
- B. Adult Occlusal/Adult Buccal**.....32803
2¹³/₁₆" (71 mm) x 7¹/₂" (190 mm) x 2⁷/₁₆" (52 mm)
- C. Adult Occlusal/Child Occlusal**.....32802
2⁷/₈" (73 mm) x 7¹/₄" (184 mm) x 2⁵/₈" (66 mm)
- D. Child Occlusal/Child Buccal**.....32806
2¹/₂" (64 mm) x 7" (178 mm) x 1¹³/₁₆" (46 mm)
- E. Lingual Upper/Lingual Lower**.....32804
1¹³/₁₆" (46 mm) x 7¹/₂" (190 mm) x 1¹³/₁₆" (46 mm)
- F. Adult Buccal/Child Buccal**.....32805
1¹³/₁₆" (46 mm) x 7³/₄" (197 mm) x 2⁷/₁₆" (62 mm)
- G. Buccal**.....32801
1¹¹/₁₆" (43 mm) x 6¹³/₁₆" (173 mm) x 1¹¹/₁₆" (43 mm)



PixSure Perfect[®]

Clinical Photography Supplies

PixSure Perfect[®]

Retract EEZ[®] Photo Lip Retractor

Innovative Design

Provides Maximum Exposure

- Features easy grip handles
- Great anatomical fit
- One size fits most
- Excellent for occlusal photographs
- Angled bend helps retract lips
- Smooth rounded edge for patient comfort
- Made in the USA

Item #: 84516 (Clear)

1 per pack

Autoclavable up to 290°F/143°C

Cold Sterilizable

Heat sterilization may cause hazing

Not recommended for Dentronix



PixSure Perfect[®]

Stainless Steel Double Ended Photo Cheek Retractors

Two Sizes in One for a Sure Fit and Full Unobstructed View

- Double ended with a large retractor on one side and a small retractor on the other
- High quality stainless steel will not discolor, yellow or become brittle like plastic
- Brushed stainless steel construction
- Accepts any method of sterilization
- Ultra-slim design for maximum comfort

Item #: DT-159

Sterilizable up to 320°F/160°C

1 pair per pack





PixSure Perfect®

Clinical Photography Supplies

PixSure Perfect®

Photo Cheek Retractors Double Ended High Density Plastic Withstands Repeated Sterilization

- Sized to comfortably fit most adults and children
- Durable plastic design
- Will not dry out, crack, or turn cloudy with proper sterilization
- Lightweight construction
- Suitable for most types of posterior photography



Item #: DT-155 (Clear)
Autoclavable up to 250°F/121°C
1 pair per pack

PixSure Perfect®

Photo Cheek Retractors Long Handles for Easy Grasp



- Durable plastic design
- Adult and child sizes
- Smooth edges for patient comfort
- Easily withstands repeated sterilization
- Provides full exposure of the posterior

A. Small (Blue)

Item #: WW-605451

B. Large (Blue)

Item #: WW-605450

Autoclavable up to 250°F/121°C

1 pair per pack



© 2021 Ortho Technology, Inc. PixSure Perfect is a registered trademark of Ortho Technology.